


BIOSAN,
Injenernaya Str., 28
630090, Novosibirsk, Russia
tel.: +7(383) 363-51-91, 363-22-40
fax: +7(383) 363-51-91
mail@biosan-nsk.ru
www.biosan-nsk.ru



DATA QUALITY SHEET

100 bp DNA Ladder (500x2) no stain



Features: The 500 bp fragment is present at increased intensity to allow easy identification

Description 100 bp DNA Ladder is ideal for determining the size of double-stranded DNA from 100 to 1,000 base pairs. The ladder consists of ten fragments that range in size from 100 to 1,000 bp in 100 bp increments. The 400 bp fragment is present at increased intensity to allow easy identification. All fragments are blunt-ended.

Fragment sizes: 1000, 900, 800, 700, 600, 500 ×2, 400, 300, 200, 100

Concentration 0.5 mg/ml

Storage 10 mM Tris-HCl (pH 8.0), 1 mM EDTA
-20°C.