


BIOSAN,
Injenernaya Str., 28
630090, Novosibirsk, Russia
tel.: +7(383) 363-51-91, 363-22-40
fax: +7(383) 363-51-91
mail@biosan-nsk.ru
www.biosan-nsk.ru



DATA QUALITY SHEET

pBlueSK/MspI DNA Marker



Features: To determine dsDNA fragment lengths ranging from 24 to 710 bp.

Description pBlueSK/MspI Marker is ideal for determining the size of double-stranded DNA from 20 to 700 base pairs. The ladder consists of 13 fragments that range in size from 24 to 710 bp. All fragments are blunt-ended.

Fragment sizes: 710, 489, 404, 328, 242, 190, 157, 147, 110, 67, 57, 34, 24.

Concentration 0.5 mg/ml

Storage 10 mM Tris-HCl (pH 8.0), 1 mM EDTA
-20°C.