


BIOSAN,  
Injenernaya Str., 28  
630090, Novosibirsk, Russia  
tel.: +7(383) 363-51-91, 363-22-40  
fax: +7(383) 363-51-91  
[mail@biosan-nsk.ru](mailto:mail@biosan-nsk.ru)  
[www.biosan-nsk.ru](http://www.biosan-nsk.ru)



## DATA QUALITY SHEET

### pBlueSK/MspI DNA Marker



**Features:** To determine dsDNA fragment lengths ranging from 24 to 710 bp.

**Description** pBlueSK/MspI Marker is ideal for determining the size of double-stranded DNA from 20 to 700 base pairs. The ladder consists of 13 fragments that range in size from 24 to 710 bp. All fragments are blunt-ended.

Fragment sizes: 710, 489, 404, 328, 242, 190, 157, 147, 110, 67, 57, 34, 24.

**Concentration** 0.5 mg/ml

**Storage** 10 mM Tris-HCl (pH 8.0), 1 mM EDTA  
-20°C.