


BIOSAN,
Injenernaya Str., 28
630090, Novosibirsk, Russia
tel.: +7(383) 363-51-91, 363-22-40
fax: +7(383) 363-51-91
mail@biosan-nsk.ru
www.biosan-nsk.ru



DATA QUALITY SHEET

DNA Marker 250 bp no stain



Features: ideal for determining the size of double-stranded DNA from 250 to 3000 base pairs

Description 250 bp DNA Ladder is ideal for determining the size of double-stranded DNA from 250 to 3000 base pairs. The ladder consists of 8 DNA fragments. All fragments are blunt-ended.

Concentration 0.5 mg/ml

Storage 10 mM Tris-HCl (pH 8.0), 1 mM EDTA
-20°C.