


BIOSAN,  
 Injenernaya Str., 28  
 630090, Novosibirsk, Russia  
 tel.: +7(383) 363-51-91, 363-22-40  
 fax: +7(383) 363-51-91  
[mail@biosan-nsk.ru](mailto:mail@biosan-nsk.ru)  
[www.biosan-nsk.ru](http://www.biosan-nsk.ru)



## DATA QUALITY SHEET

### DNA Marker 50 bp no stain



**Features:** ideal for determining the size of double-stranded DNA from 50 to 1000 base pairs

**Description** **50 bp plus DNA Ladder** is ideal for determining the size of double-stranded DNA from 50 to 1000 base pairs. The ladder consists of 23 DNA fragments. The 50, 500, 750, 1000 bp fragments are present at increased intensity to allow easy identification. All fragments are blunt-ended.

**Concentration** 0.5 mg/ml

**Storage** 10 mM Tris-HCl (pH 8.0), 1 mM EDTA  
-20°C.